

Pressure Ulcer Categorisation^{1,2}

Where Flaminal® Can Be Used


Flaminal® Hydro or Forte can be used where the tube images appear against the category. Please refer to the instructions for use on guidance on how & when to apply Flaminal®



Pack Size	PIP	NHS	Pack Size	PIP	NHS
5 x 15g tubes	324-2971	ELG021	5 x 15g tubes	324-2963	ELG022
1 x 50g tube	344-9600	ELG025	1 x 50g tube	344-9592	ELG023
500g tub	-	ELG029	500g tub	-	ELG028

Category 1: Non-Blanchable Erythema


- Intact skin with non-blanchable redness, usually over a bony prominence.
- Darker skin tones may not have visible blanching; its colour may differ from the surrounding area.
- Area may be painful, firm or soft, warmer or cooler compared to adjacent tissue.



Category 2: Partial Thickness Skin Loss


- Partial thickness loss of epidermis/dermis presenting as a shiny or dry shallow ulcer with a red/pink wound bed and without slough or bruising*
- May also present as an intact or open/ruptured serum-filled blister

* This category should not be used to describe skin tears, tape burns, perineal dermatitis, maceration or excoriation.




Category 3: Full Thickness Skin Loss

- Full thickness skin loss. Subcutaneous fat may be visible, but bone, tendon or muscle not exposed or palpable
- Slough may be present but does not obscure the depth of tissue loss. May include undermining and tunnelling.
- Depth varies by anatomical location- if no subcutaneous tissue can be shallow, but areas of significant adiposity can be extremely deep



Category 4: Full thickness Tissue Loss

- Full thickness tissue loss with exposed or palpable bone, tendon or muscle.
- Slough or eschar may be present but does not obscure depth of tissue loss. Often undermining and tunnelling.
- Depth varies by anatomical location- if no subcutaneous tissue can be shallow.
- Can extend into muscle and/or supporting structures (e.g. tendon)



Unstageable: Depth Unknown

- Full thickness tissue loss - base of the ulcer covered by slough and/or necrotic tissue
- Until enough slough and/or necrotic tissue is removed to expose the base of wound, the true depth, and category, cannot be determined.
- Dry, intact eschar on the heels should not be removed, consider referral to podiatry for advice
- Once devitalised tissue removed & category established this should be documented




Suspected Deep Tissue Injury: Depth Unknown

- Epidermis is intact, purple or maroon area of discoloured intact skin or blood-filled blister
- Area may be painful, firm, mushy, boggy, warmer or cooler compared to adjacent tissue
- May be difficult to detect in individuals with dark skin tones- may include a thin blister over a dark wound bed.
- May evolve rapidly to become covered by thin eschar.




Mucosal Pressure Ulcer:

- Develop on mucosal membranes eg. tongue, mouth, nasal passages, genital, rectum
- Cannot be categorised as the tissue does not have the same layers as the skin
- These PU are therefore unstageable (NOT unstageable)
- They are usually caused by devices and therefore should be recorded as PU (d) or mucosal pressure ulcer



Device-Related Pressure Ulcer:

- Result from the use of devices designed and applied for diagnostic/therapeutic purposes
- Some may be allocated a category of damage eg. cheeks, ears. Others may not as they appear in places that do not have the same structures as the skin eg. the mucosal membrane (see mucosal pressure ulcer)
- Where possible, a device-related ulcer should be categorised and presence of device noted by adding (d) after the category.



NOTE: The FIRST time a patient is treated for wounds which might reach the level of bones and joints or with exposed bones and joints, the patient should be under observation for at least 30 minutes after the administration of Flaminal® (Hydro or Forte)

References:

- NHS Improvement Pressure ulcer categorisation group (2019) Pressure Ulcer Categorisation. Available from <http://nhs.stopthepressure.co.uk/IMAGES TAKEN FROM THIS SOURCE>
- NHS Scotland(2021) Scottish Adaptation of the European Pressure Ulcer Advisory panel (EPUAP) Pressure Ulcer Classification Tool. Available online at: <https://www.healthcareimprovementscotland.org/idoc.ashx?docid=1f2f3549-1d38-431d-be6e-e79627e50754&version=-1>